

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

**NFR 10x**

Revision date: 25.07.2018

Product code:

Page 1 of 8

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

NFR 10x

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### Use of the substance/mixture

Laboratory use

#### Uses advised against

Any non-intended use.

### 1.3. Details of the supplier of the safety data sheet

Company name:	Dr. Lohmann Diaclean GmbH	
Street:	Oespeler Kirchweg 10	
Place:	D-44379 Dortmund	
Telephone:	+49 (0)231 - 177 285 40	
Responsible Department:	Dr. Gans-Eichler	e-mail: info@tge-consult.de
	Chemieberatung GmbH	Tel.: +49(0)251/394868-69
	Raesfeldstr. 22	www.tge-consult.de
	D-48149 Münster	

**1.4. Emergency telephone number:** +49 (0)231 - 177 285 40 Only available during office hours.

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

#### Regulation (EC) No. 1272/2008

This mixture is not classified as hazardous in accordance with Regulation (EC) No. 1272/2008.

### 2.2. Label elements

#### Additional advice on labelling

Labelling according to Regulation (EC) No. 1272/2008 [CLP]: none

### 2.3. Other hazards

The substances in the mixture do not meet the PBT/vPvB criteria according to REACH, annex XIII.  
No risks worthy of mention. Please observe the information on the safety data sheet at all times.

## SECTION 3: Composition/information on ingredients

### 3.2. Mixtures

#### Chemical characterization

The product does not contain dangerous substances according to REGULATION (EU) No. 2015/830, Annex II, Part A, 3.2.2. that must be mentioned in Chapter 3.

#### Further Information

Product does not contain listed SVHC substances > 0,1 % according to Regulation (EC) No. 1907/2006 Article 59 (REACH).

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

#### General information

In case of accident or unwellness, seek medical advice immediately (show directions for use or safety data sheet if possible).

#### After inhalation

In case of accident by inhalation: remove casualty to fresh air and keep at rest. In case of respiratory tract

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

**NFR 10x**

Revision date: 25.07.2018

Product code:

Page 2 of 8

irritation, consult a physician.

**After contact with skin**

Gently wash with plenty of soap and water. In case of skin irritation, consult a physician.

**After contact with eyes**

Rinse cautiously with water for several minutes. In case of troubles or persistent symptoms, consult an ophthalmologist.

**After ingestion**

Rinse mouth thoroughly with water. Let water be drunk in little sips (dilution effect). Do NOT induce vomiting. In all cases of doubt, or when symptoms persist, seek medical advice.

**4.2. Most important symptoms and effects, both acute and delayed**

No information available.

**4.3. Indication of any immediate medical attention and special treatment needed**

Treat symptomatically.

**SECTION 5: Firefighting measures****5.1. Extinguishing media****Suitable extinguishing media**

Carbon dioxide (CO<sub>2</sub>). Dry extinguishing powder. alcohol resistant foam. Atomized water.

**Unsuitable extinguishing media**

High power water jet.

**5.2. Special hazards arising from the substance or mixture**

Can be released in case of fire: Carbon monoxide. Carbon dioxide (CO<sub>2</sub>).Sulfur oxides. Nitrogen oxides (NO<sub>x</sub>). Gas/vapours, harmful.

**5.3. Advice for firefighters**

In case of fire: Wear self-contained breathing apparatus.

**Additional information**

Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water. Co-ordinate fire-fighting measures to the fire surroundings.

**SECTION 6: Accidental release measures****6.1. Personal precautions, protective equipment and emergency procedures**

See protective measures under point 7 and 8.

**6.2. Environmental precautions**

Discharge into the environment must be avoided.

**6.3. Methods and material for containment and cleaning up**

Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal. Clean contaminated objects and areas thoroughly observing environmental regulations.

**6.4. Reference to other sections**

Safe handling: see section 7  
Personal protection equipment: see section 8  
Disposal: see section 13

**SECTION 7: Handling and storage****7.1. Precautions for safe handling****Advice on safe handling**

Wear suitable protective clothing. (See section 8.)

**Advice on protection against fire and explosion**

Usual measures for fire prevention.

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

**NFR 10x**

Revision date: 25.07.2018

Product code:

Page 3 of 8

**Further information on handling**

General protection and hygiene measures: refer to chapter 8

**7.2. Conditions for safe storage, including any incompatibilities****Requirements for storage rooms and vessels**

Keep container tightly closed in a cool, well-ventilated place.

**Advice on storage compatibility**

Do not store together with: Explosives. Oxidizing solids. Oxidizing liquids. Radioactive substances. Infectious substances. Food and animal feedingstuff

**Further information on storage conditions**

Keep the packing dry and well sealed to prevent contamination and absorption of humidity.

Recommended storage temperature: 20°C

Protect against: Light. UV-radiation/sunlight. heat. moisture.

**7.3. Specific end use(s)**

See section 1.

**SECTION 8: Exposure controls/personal protection****8.1. Control parameters****Additional advice on limit values**

To date, no national critical limit values exist.

**8.2. Exposure controls****Appropriate engineering controls**

No special measures are necessary.

**Protective and hygiene measures**

Always close containers tightly after the removal of product. When using do not eat, drink, smoke, sniff. Wash hands before breaks and after work.

**Eye/face protection**

Wear safety glasses; chemical goggles (if splashing is possible). DIN EN 166

**Hand protection**

In case of prolonged or frequently repeated skin contact:

Wear suitable gloves.

Suitable material:

FKM (fluororubber). - Thickness of glove material: 0,4 mm

Breakthrough time  $\geq$  8 h

Butyl rubber. - Thickness of glove material: 0,5 mm

Breakthrough time  $\geq$  8 h

CR (polychloroprenes, Chloroprene rubber). - Thickness of glove material: 0,5 mm

Breakthrough time  $\geq$  8 h

NBR (Nitrile rubber). - Thickness of glove material: 0,35 mm

Breakthrough time  $\geq$  8 h

PVC (Polyvinyl chloride). - Thickness of glove material: 0,5 mm

Breakthrough time  $\geq$  8 h

The selected protective gloves have to satisfy the specifications of EU Directive 89/686/EEC and the standard EN 374 derived from it.

Before using check leak tightness / impermeability. In the case of wanting to use the gloves again, clean them before taking off and air them well.

**Skin protection**

Suitable protective clothing: Lab apron.

Minimum standard for preventive measures while handling with working materials are specified in the TRGS 500.

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**NFR 10x**

Revision date: 25.07.2018

Product code:

Page 4 of 8

**Respiratory protection**

With correct and proper use, and under normal conditions, breathing protection is not required.

**Environmental exposure controls**

No special precautionary measures are necessary.

**SECTION 9: Physical and chemical properties****9.1. Information on basic physical and chemical properties**

Physical state:	liquid	
Colour:	not determined	
Odour:	characteristic	
pH-Value:		not determined

**Changes in the physical state**

Melting point:		not determined
Initial boiling point and boiling range:		not determined
Sublimation point:		not determined
Softening point:		not determined
Pour point:		not determined
Flash point:		not determined
Sustaining combustion:		Not sustaining combustion

**Explosive properties**

none

Lower explosion limits:		not determined
Upper explosion limits:		not determined
Ignition temperature:		not determined

**Auto-ignition temperature**

Gas:

not determined

Decomposition temperature:		not determined
----------------------------	--	----------------

**Oxidizing properties**

none

Vapour pressure:		not determined
Density:		not determined
Water solubility:		not determined

**Solubility in other solvents**

not determined

Partition coefficient:		not determined
Viscosity / dynamic:		not determined
Viscosity / kinematic:		not determined
Flow time:		not determined
Vapour density:		not determined
Evaporation rate:		not determined
Solvent separation test:		not determined
Solvent content:		not determined

**9.2. Other information**

Solid content:		not determined
----------------	--	----------------

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

Revision date: 25.07.2018

**NFR 10x**

Product code:

Page 5 of 8

## SECTION 10: Stability and reactivity

### **10.1. Reactivity**

No information available.

### **10.2. Chemical stability**

The product is chemically stable under recommended conditions of storage, use and temperature.

### **10.3. Possibility of hazardous reactions**

No information available.

### **10.4. Conditions to avoid**

Protect against: UV-radiation/sunlight. heat.

### **10.5. Incompatible materials**

Materials to avoid: Oxidizing agents, strong. Reducing agents, strong.

### **10.6. Hazardous decomposition products**

Can be released in case of fire: Carbon monoxide. Carbon dioxide (CO<sub>2</sub>). Sulfur oxides. Nitrogen oxides (NO<sub>x</sub>).  
Gas/vapours, harmful.

## SECTION 11: Toxicological information

### **11.1. Information on toxicological effects**

#### **Toxicokinetics, metabolism and distribution**

No data available.

#### **Acute toxicity**

Based on available data, the classification criteria are not met.

#### **Irritation and corrosivity**

Based on available data, the classification criteria are not met.

#### **Sensitising effects**

Based on available data, the classification criteria are not met.

#### **Carcinogenic/mutagenic/toxic effects for reproduction**

Based on available data, the classification criteria are not met.

#### **STOT-single exposure**

Based on available data, the classification criteria are not met.

#### **STOT-repeated exposure**

Based on available data, the classification criteria are not met.

#### **Aspiration hazard**

Based on available data, the classification criteria are not met.

#### **Specific effects in experiment on an animal**

No data available.

## SECTION 12: Ecological information

### **12.1. Toxicity**

The product has not been tested.

### **12.2. Persistence and degradability**

The product has not been tested.

### **12.3. Bioaccumulative potential**

No indication of bioaccumulation potential.

### **12.4. Mobility in soil**

No data available.

### **12.5. Results of PBT and vPvB assessment**

## Safety Data Sheet

according to Regulation (EC) No 1907/2006

**NFR 10x**

Revision date: 25.07.2018

Product code:

Page 6 of 8

The substances in the mixture do not meet the PBT/vPvB criteria according to REACH, annex XIII.

### **12.6. Other adverse effects**

No data available.

### **Further information**

Do not allow to enter into surface water or drains.

## SECTION 13: Disposal considerations

### **13.1. Waste treatment methods**

#### **Advice on disposal**

Observe in addition any national regulations! Consult the local waste disposal expert about waste disposal.  
 Non-contaminated packages may be recycled.  
 According to EAKV, allocation of waste identity numbers/waste descriptions must be carried out in a specific way for every industry and process.  
 Control report for waste code/ waste marking according to EAKV:

#### **Waste disposal number of waste from residues/unused products**

160306 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; off-specification batches and unused products; organic wastes other than those mentioned in 16 03 05

#### **Waste disposal number of used product**

160306 WASTES NOT OTHERWISE SPECIFIED IN THE LIST; off-specification batches and unused products; organic wastes other than those mentioned in 16 03 05

#### **Waste disposal number of contaminated packaging**

150106 WASTE PACKAGING; ABSORBENTS, WIPING CLOTHS, FILTER MATERIALS AND PROTECTIVE CLOTHING NOT OTHERWISE SPECIFIED; packaging (including separately collected municipal packaging waste); mixed packaging

#### **Contaminated packaging**

Handle contaminated packages in the same way as the substance itself.

## SECTION 14: Transport information

### **Land transport (ADR/RID)**

<b><u>14.1. UN number:</u></b>	No dangerous good in sense of this transport regulation.
<b><u>14.2. UN proper shipping name:</u></b>	No dangerous good in sense of this transport regulation.
<b><u>14.3. Transport hazard class(es):</u></b>	No dangerous good in sense of this transport regulation.
<b><u>14.4. Packing group:</u></b>	No dangerous good in sense of this transport regulation.

### **Inland waterways transport (ADN)**

<b><u>14.1. UN number:</u></b>	No dangerous good in sense of this transport regulation.
<b><u>14.2. UN proper shipping name:</u></b>	No dangerous good in sense of this transport regulation.
<b><u>14.3. Transport hazard class(es):</u></b>	No dangerous good in sense of this transport regulation.
<b><u>14.4. Packing group:</u></b>	No dangerous good in sense of this transport regulation.

### **Marine transport (IMDG)**

<b><u>14.1. UN number:</u></b>	No dangerous good in sense of this transport regulation.
<b><u>14.2. UN proper shipping name:</u></b>	No dangerous good in sense of this transport regulation.
<b><u>14.3. Transport hazard class(es):</u></b>	No dangerous good in sense of this transport regulation.
<b><u>14.4. Packing group:</u></b>	No dangerous good in sense of this transport regulation.

### **Air transport (ICAO-TI/IATA-DGR)**

<b><u>14.1. UN number:</u></b>	No dangerous good in sense of this transport regulation.
<b><u>14.2. UN proper shipping name:</u></b>	No dangerous good in sense of this transport regulation.
<b><u>14.3. Transport hazard class(es):</u></b>	No dangerous good in sense of this transport regulation.

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**NFR 10x**

Revision date: 25.07.2018

Product code:

Page 7 of 8

**14.4. Packing group:** No dangerous good in sense of this transport regulation.**14.5. Environmental hazards**

ENVIRONMENTALLY HAZARDOUS: no

**14.6. Special precautions for user**

refer to chapter 6-8

**14.7. Transport in bulk according to Annex II of Marpol and the IBC Code**

not relevant

**SECTION 15: Regulatory information****15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****EU regulatory information**

2010/75/EU (VOC): No information available.  
 2004/42/EC (VOC): No information available.  
 Information according to 2012/18/EU (SEVESO III): Not subject to 2012/18/EU (SEVESO III)

**Additional information**

The mixture is classified as not hazardous according to regulation (EC) No 1272/2008 [CLP].  
 REACH 1907/2006 Appendix XVII, No (mixture): not relevant

**National regulatory information**

Water contaminating class (D): -- not water contaminating

**15.2. Chemical safety assessment**

For the following substances of this mixture a chemical safety assessment has been carried out:

**SECTION 16: Other information****Changes**

Rev. 1.0; Initial release: 25.07.2018

**Abbreviations and acronyms**

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
 CAS Chemical Abstracts Service  
 DNEL: Derived No Effect Level  
 IARC: INTERNATIONAL AGENCY FOR RESEARCH ON CANCER  
 IMDG: International Maritime Code for Dangerous Goods  
 IATA: International Air Transport Association  
 IATA-DGR: Dangerous Goods Regulations by the "International Air Transport Association" (IATA)  
 ICAO: International Civil Aviation Organization  
 ICAO-TI: Technical Instructions by the "International Civil Aviation Organization" (ICAO)  
 GHS: Globally Harmonized System of Classification and Labelling of Chemicals  
 GefStoffV: Gefahrstoffverordnung (Ordinance on Hazardous Substances, Germany)  
 LOAEL: Lowest observed adverse effect level  
 LOAEC: Lowest observed adverse effect concentration  
 LC50: Lethal concentration, 50 percent  
 LD50: Lethal dose, 50 percent  
 NOAEL: No observed adverse effect level  
 NOAEC: No observed adverse effect level  
 NTP: National Toxicology Program  
 N/A: not applicable  
 OSHA: Occupational Safety and Health Administration  
 PNEC: predicted no effect concentration  
 PBT: Persistent bioaccumulative toxic

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**NFR 10x**

Revision date: 25.07.2018

Product code:

Page 8 of 8

RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail )

SARA: Superfund Amendments and Reauthorization Act

SVHC: substance of very high concern

TRGS Technische Regeln fuerGefahrstoffe

TSCA: Toxic Substances Control Act

VOC: Volatile Organic Compounds

VwVwS: Verwaltungsvorschrift wassergefaehrdender Stoffe

WGK: Wassergefaehrdungsklasse

**Further Information**

Classification according EC regulation 1272/2008 (CLP): - Classification procedure:

Health hazards: Calculation method.

Environmental hazards: Calculation method.

Physical hazards: On basis of test data and / or calculated and / or estimated.

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

---

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*